

Disinfectants, B/Braun

① HELIPUR® H plus N

Formaldehyde-free concentrate for disinfecting thermostable and nonthermo-stable materials. Evaluated using the new DGHM (German Hygiene and Microbiology Society) method. Formaldehyde-free. Suitable for surgical instruments and non-thermostable items including flexible endoscopes. Can be used in an ultrasonic bath. Good material compatibility. Effective across a wide spectrum (bacteria incl. TB, fungus, viruses incl. HBV and HIV). Concentrations for use acc. to DGHM: 1% = 60 min, 1% = 30 min, 1.5% = 15 min. Biodegradable.

HELIPUR® H plus N Cleaner for severe soiling and incrustations.

Article	Vol.	Order No.
Dispensing bottle	1000 ml	35 42 31110
Canister	5 l	35 42 31510
HELIPUR H plus N Cleaner dispenser bottle	1000 ml	35 42 31610

② MELISEPTOL New Formula®

Ready-to-use alcoholic disinfectant for spraying or wiping. For medical inventory and other alcohol-resistant surfaces. Contains 44 g ethanol and not other active ingredients, is free of quaternary ammonium compounds, amines, and aldehydes. Effective against bacteria (incl. MRSA, TbB and mycobacteria) and fungus; limited-virucidal (incl. HBV, HCV, HIV) and effective against rotaviruses, adenoviruses and noroviruses. DGHM-/VAH-listed and in the IHO virucide list. Impact time: 1 minute.

Article	Vol.	Order No.
Spray bottle (with spray head)	250 ml	10 42 19763
Spray bottle (without spray head)	1000 ml	10 42 19764
Dosage bottle	1000 ml	10 42 19761
Canister	5 l	10 42 19758
Spray head, red		35 42 32900

③ MELISEPTOL® Foam pure

Fast-acting foam for disinfection of small surfaces. Sprays without formation of aerosol, free of perfume, active against bacteria incl. TbB, limited action against viruses. Well tolerated also by material susceptible against alcohol, e. g. medical devices, acrylic glass, examination couches.

Article	Volume	Order No.
Canister	5 l	10 42 19288
Spray bottle (with spray head)	750 ml	10 42 19171

④ MELISEPTOL® rapid

Quick-acting, ready-to-use disinfectant for your medical inventory. MELISEPTOL® rapid has a reaction time of just 1 minute against bacteria incl. tuberculosis pathogens, fungi incl. moulds, as well as against HBV, HCV, HIV, rota, adeno and vaccina viruses.

Article	Vol.	Order No.
Spray bottle (with spray head)	250 ml	35 42 32600
Dispensing bottle (without spray head)	250 ml	35 42 32700
Hand spray bottle	1000 ml	35 42 32800
Canister	5 l	35 42 32805

⑤ MELISEPTOL® HBV wipes

Sealed in aluminium foil to prevent drying out. Areas of use same as MELISEPTOL®.

Article	Order No.
Dispenser box with 100 wipes	35 42 32902
Refill pack with 100 wipes	35 42 32913

⑥ HEXAQUART® forte

Aldehyde-free disinfectant concentrate with high cleaning power for inventory and floors.

Concentration: 1% / 30 min, 0.5% / 1 h, 0.25% / 4 h

Article	Vol.	Order No.
Dispensing bottle	1000 ml	10 42 19105
Canister	5 l	10 42 19106
Dispensing pump for 5-l canister, stroke 15 ml (for 1.5%)		35 42 32911
Dispensing pump for 5-l canister, stroke 20 ml (for 0.5%)		35 42 32912

⑦ MELSEPT® SF

Formaldehyde-free surface disinfectant for the disinfection and cleaning of washable surfaces and floors; bactericidal, fungicidal, tuberculocidal, virucidal. Effective against Hepatitis-B viruses.

Concentration: 0.5% = 60 min.

Article	Vol.	Order No.
Dispensing bottle	1000 ml	35 42 33100
Canister	5 l	35 42 33500

Dispensing beaker (free)

