

MBT Adhesive Indicator Tapes, IDL

For monitoring sterilization processes. Reliable and positive identification of sterilization quality.



WW7R3068



Article	Order No.
---------	-----------

- | | |
|---|--------------------|
| ① Stericlavable indicator tape, 1 reel of 55 m, 18 mm wide | 29 78 55191 |
|---|--------------------|

The indicator on the tape changes from light beige to dark brown during steam sterilization in autoclaves. It also reacts on the action of pressure and temperature after a minimum time - at 120°C after 20 minutes and at 135°C after 7 minutes.

Application: Sensitive sterilization products in steam sterilizers, e.g. culture solutions and the like.

Article	Order No.
---------	-----------

- | | |
|---|--------------------|
| ② Hot air indicator tape, 1 reel of 55 m, 18 mm wide | 29 78 55192 |
|---|--------------------|

The indicator in the tape changes from light green to brown during hot-air sterilization. It also reacts to temperature and time - at 160°C after 60 minutes and at 180°C after 30 minutes.

Applications: Non-sensitive equipment and aids.

Article	Order No.
---------	-----------

- | | |
|--|--------------------|
| ③ Hand dispenser for indicator tape | 29 78 55193 |
|--|--------------------|

Stable. For reliable dispensing and cutting of the indicator tape.

