

Distilling Columns, Rettberg

① Hempel type columns without glass jacket

for simple fractioning, with packing carrier, borosilicate glass 3.3.

Gr. joint	I.D.	Filling H	Order No.
NS	mm	mm	
14.5/23	15	200	13 40 65214
14.5/23	15	300	13 40 65314

Gr. joint	I.D.	Filling H	Order No.
NS	mm	mm	
29/32	30	300	13 40 65329
29/32	30	500	13 40 65529

② Hempel type columns with slip-on glass jacket

for simple fractioning, with packing carrier, borosilicate glass 3.3.

Gr. joint	I.D.	Filling H	Order No.
NS	mm	mm	
14.5/23	15	200	13 40 65614
14.5/23	15	300	13 40 65714

Gr. joint	I.D.	Filling H	Order No.
NS	mm	mm	
29/32	30	300	13 40 65829
29/32	30	500	13 40 65929

③ Vigreux columns with slip-on glass jacket

Borosilicate glass 3.3

Gr. joint	I.D.	Filling H	Order No.
NS	mm	mm	
14.5/23	15	200	13 40 66215
14.5/23	15	300	13 40 66315

Gr. joint	I.D.	Filling H	Order No.
NS	mm	mm	
29/32	25	300	13 40 66330
29/32	25	600	13 40 66630

