

Filter Funnels, DURAN®, DURAN Group / IDL



WWEQ2076



Filter funnels are connected to the filtering flask by means of a cone-shaped rubber seal (GUKO).

Porosity	Vol. ml	Plate Ø mm	Stem O.D. mm	Order No.
1	50	35	10	11 25 85201
2	50	35	10	11 25 85202
3	50	35	10	11 25 85203
4	50	35	10	11 25 85204
5	50	35	10	11 25 85205
1	75	45	10	11 25 85211
2	75	45	10	11 25 85212
3	75	45	10	11 25 85213
4	75	45	10	11 25 85214
5	75	45	10	11 25 85215
1	125	60	10	11 25 85221
2	125	60	10	11 25 85222
3	125	60	10	11 25 85223
4	125	60	10	11 25 85224
5	125	60	10	11 25 85225
1	500	95	22	11 25 85231
2	500	95	22	11 25 85232
3	500	95	22	11 25 85233
4	500	95	22	11 25 85234
5	500	95	22	11 25 85235
1	1000	120	22	11 25 85241
2	1000	120	22	11 25 85242
3	1000	120	22	11 25 85243
4	1000	120	22	11 25 85244
5	1000	120	22	11 25 85245
1	4000	175	30	11 25 85261
2	4000	175	30	11 25 85262
3	4000	175	30	11 25 85263
4	4000	175	30	11 25 85264

