

Apparatus for Determination of Essential Oils in Drugs, Starke

Article	Order No.
Apparatus for determining ethereal oils	14 59 50129

In accordance with DAB, distillation head with ground joint-cone size 29/32, Allihn condenser, vent plug with ground joint size 7/16, measuring tube 1 ml: 0.01 with 3 ml ball, Czako stopcock, Teflon plug, without piston

Article	Order No.
Apparatus for determining ethereal oils	14 59 50229

Same as above but Czako stop with glass plug

On request, measuring pipes are also available with different contents, ball, or measuring pipe.



STARKE
Laborbedarf
webcode  **WWFNK005**

