

Water Analysis Apparatuses, Starke

Water Analysis Apparatuses

Article	Order No.
Oxygen bottle according to Winkler	14 59 40010

To determine the oxygen dissolved in water, complete with stopper and matte-finish marking area, capacity exactly calibrated. Capacity: 100 to 150 ml

Article	Order No.
① Oxygen bottle according to Winkler	14 59 40011

Same as above. Capacity: 250 to 300 ml

Article	Order No.
② Sedimenting vessel, Imhoff type	14 59 40013

Without stopcock, cap. 1000 ml, graduated: 0 to 2 ml in 0.1 ml; 2 to 10 ml in 0.5 ml; 10 to 40 ml in 1 ml; 40 to 100 ml in 2 ml, with mark at 1000 ml

Article	Order No.
Sedimenting vessel, Imhoff type	14 59 40014

With stopcock with ground joint and solid glass plug, with threaded stopcock lock, cap. 1000 ml, graduated: 0 to 2 ml in 0.1 ml; 2 to 10 ml in 0.5 ml; 10 to 40 ml in 1 ml; 40 to 100 ml in 2 ml, with mark at 1000 ml

Article	Order No.
③ Sedimenting vessel, Imhoff type	14 59 40015

Without stopcock, cap. 1000 ml, graduated: 0 to 2 ml in 0.1 ml; 2 to 10 ml in 0.5 ml; 10 to 40 ml in 1 ml; 40 to 100 ml in 2 ml; 100 to 1000 ml in 100 ml

Article	Order No.
Sedimenting vessel, Imhoff type	14 59 40016

With stopcock with ground joint and solid glass plug, with threaded stopcock lock, cap. 1000 ml, graduated: 0 to 2 ml in 0.1 ml; 2 to 10 ml in 0.5 ml; 10 to 40 ml in 1 ml; 40 to 100 ml in 2 ml; 100 to 1000 ml in 50 ml

Article	Order No.
④ BOD flask	14 59 40018

Stopper with ground joint 18.8, capacity with stopper inserted: 100 ml

- ⑤ **Overflow graduated flask**
cut off a capacity mark, plane-ground.

Volume ml	Order No.	Volume ml	Order No.
21.7	10 59 40025	164	14 59 40164
22.7	14 59 40027	200	14 59 40200
43.5	10 59 40043	250	14 59 40250
56	14 59 40056	300	14 59 40300
94	14 59 40094	360	14 59 40360
97	14 59 40097	365	14 59 40365
100	14 59 40100	428	14 59 40428
150	14 59 40150	432	14 59 40432
157	14 59 40157		

