

SICCO Glove Box Air Stream with Transfer Chamber, Bohlender

Aluminium frame with panels made of acrylic glass, base panel made of polyethylene. Front panel with two glove ports and one pair of natural rubber gloves size 9, two lead-ins for cables on the back side (upper left and upper right). Transfer chamber with two doors on the right, external door with latch to prevent accidental opening, the tray extends into the workspace with opening of the internal door for easy removal of inserted products. Integrated HEPA filtering unit in the top side for a high-purity atmosphere in the workspace. Manometer to display the degree of filter saturation. Pipe sockets made of PE are included in delivery for converting to the different operating modes: internal air circulation, fresh air supply and exhaust. The conversion can be made without any tooling. Exchangeable filter cartridge. Maintenance-free and trouble-free motor. Electrical power connection 230 V / 50 Hz is required.

Technical data:

Panels: clear
Material: PE
Temperature resistance: from -20°C to +70°C
Use: normal pressure
Stand: rubber feet
Maximum all-over load Box: 40 kg
Usable volume Box: 290 l

Outer dimensions W x H x D	Wght	Usable interior space W x H x D	Usable interior transfer chamber W x H x D	Hinged lid opening W x H	Order No.
mm	kg	mm	mm	mm	
1210x1015x660	42	850x650x540	220x275x320	850x425	10 17 64508

