

Oak Ridge Centrifuge Tubes with Screw Sealing Cap Assemblies, NALGENE®, IDL

Centrifuge tubes with leak-proof, autoclavable screw closures made of polypropylene, white and integrated seals. Closures size 16-ml and bigger with handle.

① Centrifuge tubes made of polycarbonate

Volume ml	Order No.	PU pc/Pack
10	23 87 31381	10
16	23 87 31380	10
30	23 87 31383	10
50	23 87 31384	10

❶ Available only in original packaging units. No fractions possible.

② Centrifuge tubes made of polypropylene-copolymer

Volume ml	Order No.	PU pc/Pack
10	23 87 31391	10
16	23 87 31390	10
30	23 87 31393	10
50	23 87 31394	10

❷ Available only in original packaging units. No fractions possible.

③ Centrifuge tubes made of polysulfone

Polypropylene screw cap; silicone seal

Same stability and sturdiness as polycarbonate but much better resistance to acids and alkali. Withstand up to 50,000 × g in refrigerated or non-refrigerated centrifuges. With sealing cap closure option. Unscrew the cap before autoclaving and place it on the neck of the tube. The tube must always be no less than 80% full.

Volume ml	Order No.	PU pc/Pack
50	10 87 31374	10

❸ Available only in original packaging units. No fractions possible.

