

Carboy Accessories, Thermo Scientific Nalgene, IDL



WWRJN

thermo scientific

① Filling and venting closures

Non-sterile, made of polypropylene with TPE seal, TPE caps as well as 2x Nalgene® 550 platinum-hardened silicone tubes (inside and outside).

Cap Ø mm	Tubing connectors	Outer Ø mm	Order No.
53	3	6	22 87 21618
83	3	6	22 87 21619
83	3	12	22 87 21620

Tube fitting

Article	Order No.	PU pc/Pack
---------	-----------	---------------

② Tube fitting with tube connector and thread

22 87 21621

2

made of polypropylene with silicone intermediate seal and acetate union nut for screw connections or plastic containers for subsequent installation (16 mm bore required). Screwed cap not comprised in the delivery of the retrofit kit.

❗ Available only in original packaging units. No fractions possible.

Quick-release stopcocks

Article	Order No.
---------	-----------

③ Quick-release stopcock made of Tefzel ETFE

22 87 21622

for tubes with 16 mm inner Ø, suitable for all Nalgene® canisters with threaded connection.

Article	Order No.
---------	-----------

④ Quick-release stopcock made of polypropylene

22 87 21623

with PTFE seal and plug-on tube connector with 6 mm outer Ø, suitable for all Nalgene® canisters with threaded connection.

Spigot Closure

Spigot closure for all Nalgene balloon bottles with cock, made of white polypropylene with TPE plug; enclosed with delivery of all Nalgene balloon bottles with drain cock.

Article	Order No.	PU pc/Pack
Spigot closure	10 87 30010	4

