

Distilling Apparatus, Rettberg



WWQVD136

rettberg.biz
GLASAPPARATEBAU • LABORATORY GLASSWARE

Distilling apparatus, complete

Made of borosilicate glass 3.3, consisting of the components listed below.

Article	Order No.
Distilling apparatus, complete unit	13 70 75000

Individual components

Article	Order No.
Distilling link/Claisen type	13 40 75014
with fused Liebig condenser, 2 standard ground sockets 14.5, 2 × standard ground cones 29, jacket length approx. 330 mm (fitted length approx. 250 mm), condensing up to below the cone at the outlet, condensing water connections with unscrewable plastic hose connections (PP)	
Intermediate receiver	13 40 91424
Anschütz-Thiele type, with inclined stub, capacity 50 ml, graduated, standard ground cone and socket 29, vacuum connection with unscrewable plastic hose connections (PP).	
Round bottom flask 250 ml, NS 29	13 40 20236
Evaporator flask	13 40 24254
1 ltr., middle neck ground joint 29, side neck ground joint 14.5, inclined	
Stopper with conical ground joint NS 14.5	13 40 39114
Gas supply tube, NS 14.5 with PTFE valve	13 40 37211
Thermometer with standard ground joint	10 70 00128
with standard ground joint 14.5, 0 to +250°C, fitted length 75 mm	
Ground retaining device, NS 29 (3 pieces needed)	13 70 12029
Heater 1000 ml	13 70 21110
Scissor-type table	13 70 75060
Stand set	13 70 75070

