

Heating Plates behrotest® with Metal Safety Grid, Behr

Heating Plates behrotest® with Metal Safety Grid

Adjustable hotplates for general laboratory use. A metal grid protects the user from inadvertent contact with the hot surface.

Article	Ø mm	Heating cap. W	Mains connec.	Order No.
① Heating plate behrotest® KP 1	94	500	230 V	10 48 00504
Heating plate behrotest® KP 2	145	1100	230 V	10 48 00514
Heating plate behrotest® KP 3	220	2000	230 V	10 48 00511
Heating plate behrotest® KP 4	180	2000	230 V	10 48 00513

Top molds for KP 1

Article	Order No.
② Top mold AM 100/94 for 100 ml round bottom flask incl. distance inlay	10 48 17701
Top mold AM 250/SET for 250 ml round bottom flask with distance stays	10 48 17706
Top mold AM 500/SET for 500 ml round bottom flask with distance stays	10 48 17707

Top molds for KP 2

Article	Order No.
③ Top mold AM 1000/145 for 1000 ml round bottom flask incl. distance inlay	10 48 17705

Serial heaters behrotest®

Serial heaters. Each heater is controlled individually and has a control LED. With holders for the telescope bars available as an option.

Article	Order No.
④ Serial heater HB 4 4 individually adjustable units, Ø 94 mm	10 48 00506
Serial heater HB 6 6 individually adjustable units, Ø 94 mm	10 48 00507
Serial heater HB 8 8 individually adjustable units, Ø 94 mm	10 48 00512
Holder for HB 4 incl. 4 support rods	10 48 00508
Holder for HB 6 incl. 6 support rods	10 48 00509
Holder for HB 8 incl. 8 support rods	10 48 00510

Top molds for serial heaters HB 4, HB 6 and HB 8

Maximum capacity in serial heater HB 4: 4 top molds
Maximum capacity in serial heater HB 6: 6 top molds
Maximum capacity in serial heater HB 8: 8 top molds

Article	Order No.
② Top mold AM 100/94 for 100 ml round bottom flask incl. distance inlay	10 48 17701
Top mold AM 250/SET for 250 ml round bottom flask with distance stays	10 48 17706
Top mold AM 500/SET for 500 ml round bottom flask with distance stays	10 48 17707



①



②



③



④

