

schuett petriturn Turntables for Petri Dishes, schuett-biotec

For uniform plating and inoculating of liquids on agar surfaces in petri dishes. Easy to handle, antifatique. Made in Germany.

① schuett petriturn-M

Manual operation, constant rotary motion. Powerful and persistent due to high mass (weight) of turntable. Usable on both sides, that means two work spaces in one device by easiest turn-over. For petri dishes with dia. 90 up to 150 mm.

Article	Ø × H mm	Order No.
schuett petriturn-M, manual	160 × 45	88 31 31001

② schuett petriturn-E

Electrically driven. Constant rotary motion. Simple contactless handling of the turntable, started via IR sensor (operating time adjustable from 3 to 120 sec), permanent running, connection of a foot-switch (optional). Continuously adjustable speed of 10 to 120 rpm. For petri dishes up to dia. 100 mm (dia. 150 mm optional). Compact stainless steel casing.

Article	Ø × H mm	Order No.
schuett petriturn-E, electrical, 230 V	160 × 70	88 31 31011
Foot switch for schuett petriturn-E		88 31 31013
Turntable-attachment for petri-dishes to Ø 150 mm		10 31 81015
Instrument tray, stainless steel		10 31 88352

Turntables for petri dishes of other diameters available on request.

Inoculating hooks for plating (Drigalski spatula)

Article	L mm	W mm	Order No.
③ Inoculating hook, glass	145	40	88 31 31014
④ Inoculating hook, stainless steel	180	58	88 31 31012

